

For Research Use Only

# Myc-Trap® 2.0 Agarose, Kit for Immunoprecipitation



[www.ptgcn.com](http://www.ptgcn.com)

Catalog Number: yt2ak

## Basic Information

**Catalog Number:**

yt2ak

**Applications:**

IP, Co-IP

**Conjugate:**

Agarose beads; ~90 um (cross-linked 4% agarose beads)

**Host:**

Alpaca

**Type:**

Nanobody

**Class:**

Recombinant

## Description

The ChromoTek Myc-Trap® Agarose, Kit for Immunoprecipitation consists of an anti-Myc NANOBODY®/VHH, which is coupled to agarose beads. It also contains lysis, wash, and elution buffers that can be used for the immunoprecipitation of Myc-fusion proteins from cell extracts of various organisms.

## Specificity/Target

Binds specifically to the Myc-tag (sequence EQKLISEEDL) at the N-terminus, C-terminus, or internal site of the fusion protein. Endogenous c-myc is NOT bound.

## Elution buffer

2x SDS-sample buffer (Lämmli), 200 mM glycine pH 2.5, 0.1 mg/ml ChromoTek 2x Myc-peptide (2yp) in PBS pH 7.4

## Affinity (K<sub>D</sub>)

770 nM

## Storage

**Storage:**

Shipped at ambient temperature. Upon receipt store at +4°C. Stable for one year. DO not freeze!

**Storage Buffer:**

20% ethanol

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [www.ptgcn.com](http://www.ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunoprecipitation of GFP-Myc fusion protein from HEK293T cells using Myc-Trap® 2.0 Agarose. IP was done using both untransfected (mock) and transfected (GFP-myc) cells. I: Input; F: Flow-through, B: Bound.



Comparison of pulldown efficacy between the Myc-Trap® 2.0 Agarose (left) and the original Myc-Trap Agarose (right). Both products were used to immunoprecipitate GFP-myc fusion proteins from untransfected (mock) and transfected (GFP-myc) HEK293T cells. The Myc-Trap 2.0 has higher affinity for myc-tagged proteins and is able to pulldown more GFP-Myc protein than the Myc-Trap.