

For Research Use Only

IHCeasy PSME2 Ready-To-Use IHC Kit

Catalog Number: KHC0449

General Information

Sample type: FFPE tissue Cited sample type: Reactivity: Human, Mouse, Rat Cited Reactivity: Assay type: Immunohistochemistry Primary antibody type: Rabbit Polyclonal Secondary antibody type: Polymer-HRP-Goat anti-Rabbit

Kit Component

| Component | Size | Concentration |
|--------------------------|--------------------|---------------|
| Antigen Retrieval Buffer | 100 mL | 50× |
| Washing Buffer | 100 mL ×2 | 20× |
| Blocking Buffer | 5 mL | RTU |
| Primary Antibody | 5 mL | RTU |
| Secondary Antibody | 5 mL | RTU |
| Chromogen Component A | 0.2 mL | RTU |
| Chromogen Component B | 4 mL | RTU |
| Signal Enhancer | 5 mL | RTU |
| Counter Staining Reagent | 5 mL | RTU |
| Mounting Media | 5 mL | RTU |
| Control Slide | 1 slide (Optional) | FFPE |
| Datasheet | 1 Сору | |
| Manual | 1 Сору | |

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.

Background


The principal function of the proteasome is targeted degradation of intracellular proteins. Activity of the 20S proteasome is controlled by regulatory complexes that bind to the ends of the cylindrical proteasome. 11S regulator (REG or PA28), is a complex that is thought to activate proteasomes toward the production of antigenic peptides. Human PSME1 and PSME2 genes encode the two proteasome activators PA28 alpha and beta, respectively, which have been implicated in antigen processing for loading class I MHC molecules. The PA28 activator complex enhances the generation of class I binding peptides by altering the cleavage pattern of the proteasome.

Synonyms

PA28B, PA28beta, PSME2, REG beta, REGbeta


For technical support and original validation data for this product please contact: T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data




Immunohistochemical analysis of paraffinembedded human liver tissue slide using KHC0449 (PSME2 IHC Kit).

Immunohistochemical analysis of paraffinembedded human liver cancer tissue slide using KHC0449 (PSME2 IHC Kit).




Immunohistochemical analysis of paraffinembedded human kidney tissue slide using KHC0449 (PSME2 IHC Kit).



Immunohistochemical analysis of paraffinembedded human renal cell carcinoma tissue slide using KHC0449 (PSME2 IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse liver tissue slide using KHC0449 (PSME2 IHC Kit).



Immunohistochemical analysis of paraffinembedded rat liver tissue slide using KHC0449 (PSME2 IHC Kit).