



IHCeasy TAOK3 Ready-To-Use IHC Kit

Catalog Number: KHC2908

General Information

Sample type: FFPE tissue Cited sample type: Reactivity: Human, Mouse, Rat Cited Reactivity: Assay type: Immunohistochemistry Primary antibody type: Mouse Monoclonal

Secondary antibody type: Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

TAOK3 (Serine/threonine-protein kinase TAO3) is also named as DPK, JIK, KDS, MAP3K18 and belongs to the STE Ser/Thr protein kinase family. It acts as an activator of the p38/MAPK14 stress-activated MAPK cascade. In response to DNA damage, it is involved in the G2/M transition DNA damage checkpoint by activating the p38/MAPK14 stress-activated MAPK cascade, probably by mediating phosphorylation of upstream MAP2K3 and MAP2K6 kinases. And this protein can be autophosphorylated.


Synonyms

JIK, CTCL-associated antigen HD-CL-09, Cutaneous T-cell lymphoma-associated antigen HD-CL-09, Dendritic cell-derived protein kinase, DPK


Selected Validation Data




Immunohistochemical analysis of paraffinembedded human stomach cancer tissue slide using KHC2908 (TAOK3 IHC Kit).




Immunohistochemical analysis of paraffinembedded mouse spleen tissue slide using KHC2908 (TAOK3 IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse kidney tissue slide using KHC2908 (TAOK3 IHC Kit).



Immunohistochemical analysis of paraffinembedded mouse stomach tissue slide using KHC2908 (TAOK3 IHC Kit).



Immunohistochemical analysis of paraffinembedded rat kidney tissue slide using KHC2908 (TAOK3 IHC Kit).