



# PRODUCT DESCRIPTION

FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS OR ANIMALS

ProteinTech Group, Inc  
2201 W. Campbell Park Dr.  
Chicago, IL 60612  
Toll Free: 1-888-4PTGLAB  
In Illinois: 1312-455-8498  
Fax: 1-312-455-8408  
Web: www.ptglab.com


## Mouse anti- ACTIN, BETA Monoclonal Antibody

BC014861

**ACTIN, BETA**  
**ACTB,BETA-ACTIN**

### The antibody:

Format : Caprylic Acid and  
Ammonium Sulfate purified  
Catalog No. : 60008-1  
Clone No. : 7D2C10  
Weight/Volume: 100 µg/150 µl  
Tested Applications: WB, ELISA



Human 293 cell  
lysates were  
subjected to 10%  
SDS PAGE  
followed by  
western blot with  
1:1000

### The production

Immunogen: Hu, recombinant protein.  
Specificity: Human  
Host: Mouse  
Antibody species: Monoclonal

### Formulation:

100µg of Caprylic Acid and Ammonium Sulfate  
purified IgM in 150µl of 70% storage buffer (0.1M  
Tris-HCl, pH 7.4, 0.15M NaCl, 0.04% merthiolate)  
and 30% glycerol.

### Handling Recommendations:

Aliquot into microcentrifuge tubes and store at -  
20°C. Avoid repeated freeze/thaw cycles, which  
may damage IgG and affect product performance.