

For Research Use Only

CCDC83 Polyclonal antibody

Catalog Number: 20902-1-AP



Basic Information

Catalog Number:

20902-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15023

GenBank Accession Number:

BC040280

GeneID (NCBI):

220047

UNIPROT ID:

Q8IWF9

Full Name:

coiled-coil domain containing 83

Calculated MW:

413 aa, 49 kDa

Observed MW:

49 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Jurkat cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Jurkat cells were subjected to SDS PAGE followed by western blot with 20902-1-AP (CCDC83 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.