

For Research Use Only

ASTL Polyclonal antibody

Catalog Number: 21418-1-AP



Basic Information

Catalog Number:

21418-1-AP

Size:

240 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15997

GenBank Accession Number:

BC107127

GeneID (NCBI):

431705

UNIPROT ID:

Q6HA08

Full Name:

astacin-like metallo-endopeptidase
(M12 family)

Calculated MW:

431 aa, 46 kDa

Observed MW:

43 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HEK-293 cells were subjected to SDS PAGE followed by western blot with 21418-1-AP (ASTL Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.