

For Research Use Only

RyDEN Polyclonal antibody

Catalog Number: 27865-1-AP

Featured Product

4 Publications



Basic Information

Catalog Number:

27865-1-AP

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26838

GenBank Accession Number:

BC035817

GeneID (NCBI):

55337

UNIPROT ID:

Q9NUL5

Full Name:

chromosome 19 open reading frame
66

Observed MW:

38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse, Rat

Cited Species:

human, mouse, monkey

Positive Controls:

WB : THP-1 cells, mouse liver tissue, rat liver tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Matthias M Zimmer	34893599	Nat Commun	WB
Natasha W Hanners	34873063	Proc Natl Acad Sci U S A	WB
Inès Bribes	40658726	PLoS Negl Trop Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



THP-1 cells were subjected to SDS PAGE followed by western blot with 27865-1-AP (RyDEN antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 27865-1-AP (RyDEN antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.