

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Affinipure Goat Anti-Mouse IgG, Fcγ Subclass 1 Specific

Catalog Number: SA00012-1



## Information

Catalog Number:  
SA00012-1  
Size:  
50ul/100ul  
Applications:  
ELISA, WB

Product name:  
HRP-conjugated Affinipure Goat Anti-Mouse  
IgG, Fcγ Subclass 1 Specific  
Physical State:  
Liquid  
Conjugation:  
Horseradish Peroxidase

## Recommended Dilutions

1:2000-1:40,000 for western blotting with ECL substrates  
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:  
Store at -20°C. Avoid repeated freeze / thaw cycles.  
Storage Buffer:  
0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6  
Stabilizer:  
3mg/ml BSA  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.