

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Goat Anti-Mouse IgA ( $\alpha$ )

Catalog Number:SA00012-7



## Information

Catalog Number:  
SA00012-7  
Host:  
Goat  
Applications:  
ELISA, WB, Dot blot

Reactivity:  
Mouse  
Physical State:  
Liquid  
Conjugation:  
HRP

## Recommended Dilutions

1:250-1:1500 for ELISA  
1:1500-1:5,000 for Western blotting with ECL substrate

## Stabilizer

3 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

