



PRODUCT DESCRIPTION

**FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS**

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com

LMNB1 PolyClonal Antibody

Catalog No. : 12987-1-AP

Size : 37ug / 150ul

Source : Rabbit


IsoType : IgG

Synonyms :


LMNB1;

LMNB1;LMNB1;ADLD;LMN;L


MN2;LMNB;MGC111419 ;



Hela cell were subjected to SDS PAGE followed by western blot with 12987-1-AP(LMNB1 Antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded Breast cancer using 12987-1-AP(LMNB1 Antibody) at Dilution 1:100 (under 10x lens)



Immunofluorescent analysis of HepG2 cells, using LMNB1 antibody 12987-1-AP at 1:50 dilution and FITC-labeled donkey anti-rabbit IgG(green).

Applications

Tested Applications:

ELISA,WB,IHC,IF

Specificity Reaction:

Human,Mouse,Rat; Other species are not tested.Please decide the specificity by homology

Positive WB detected In:

Hela cell,Hepg2 cell;Human brain;Human kidney;Jurkat cell;K562 cell;Mouse eye

Calculated LMNB1 MW:

66kd

Observed LMNB1 MW:

66kd

Positive IHC detected In:

Breast cancer

Application key: WB=Western Blotting ,IHC=Immunohistochemistry, IF=Immunofluorescence

Immunogen Information

Immunogen: ag3631

GenbankNo: BC012295

GeneID: 4001

Full Name: lamin B1

Product Information

Purify: Antigen Affinity Purified

Storage: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.