



PRODUCT DESCRIPTION

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS

Antibody to MBP

Catalog Number: 15089-1-AP

Isotype: IgG

Source: Rabbit

Product Information

Description: Purified Rabbit Anti MBP Polyclonal antibody

Immunogen: Recombinant Protein

Clonality: Polyclonal

Purity: Antigen Affinity Purified

Form: Liquid

Storage Condition: -20°C, Avoid freeze / thaw cycles.

Storage Buffer: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.

Concentration: 400ug/ml

Applications


Tested Applications: WB; IP; ELISA

Recommended Titers: ELISA 1:1,000~100,000

Western blotting 1:1000~3,000

Specificity: Maltose binding protein tag.

Positive control: MBP-tagged protein



MBP-tagged protein were subjected to SDS PAGE followed by western blot with affinity purified Ab at dilution of 1:1000-3000