



PRODUCT DESCRIPTION

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS

Antibody to Flag Tag

Catalog Number: 60002-1-Ig

Isotype: IgG2b

Source: Mouse

Product Information

Description: Purified Mouse Anti-Flag monoclonal antibody

Immunogen: DYKDDDDK(FLAG) synthetic peptide conjugate to KLH

Clonality: Monoclonal

Purification: Protein A Affinity Purified from mouse ascites fluid

Storage Buffer: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.

Storage Condition: Store at -20°C, Avoid freeze / thaw cycles.

Concentration: 1mg/ml

Size: 100ul

Applications


Tested Applications: WB; IP; ELISA ,Dot

Recommended Titers: ELISA 1:8000~100,000

Western blotting 1:3000~8000

Specificity: Recognizes the N-terminal ,C-terminal or internal DYKDDDDK-tagged fusion proteins

Positive control: Flag-tagged protein



Flag-tagged fusion protein(32kD) were subjected to SDS-PAGE followed by Western Blot with 60002-1-Ig at dilution of 1:6000.