



# PRODUCT DESCRIPTION

ProteinTech Group, Inc  
2201 W. Campbell Park Dr.  
Chicago, IL 60612  
Toll Free: 1-888-4PTGLAB  
In Illinois: 1312-455-8498  
Fax: 1-312-455-8408  
Web: www.ptglab.com

FOR IN VITRO RESEARCH USE ONLY  
NOT FOR USE IN HUMANS OR ANIMALS

## Antibody to GFP

**Catalog Number:** 66002-1-Ig

**Isotype:** IgG2a

**Source:** Mouse

### Product Information

**Description:** Purified Mouse Anti-GFP monoclonal antibody

**Immunogen:** Recombinant Protein

**Clonality:** Monoclonal

**Purify:** Protein A Affinity Purified from mouse ascites fluid

**Form:** Liquid

**Storage Condition:** -20°C, Avoid freeze / thaw cycles.

**Storage Buffer:** PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.

**Concentration:** 0.9mg/ml

**Size:** 100ul

### Applications


**Tested Applications:** WB; IP; ELISA

**Recommended Titers:** ELISA 1:10,000~100,000

Western blotting 1:5,000~10,000

**Specificity:** Reacts with *Aequorea Victoria*.

**Positive control:** GFP-tagged protein



GST-tagged GFP fusion protein (57kd, 20ng/lane) was subjected to SDS-PAGE followed by Western Blot with 66002-1-Ig at dilution of

Lane 1 1:10000	Lane 5 1:160000
Lane 2 1:20000	Lane 6 1:320000
Lane 3 1:40000	Lane 7 1:640000
Lane 4 1:80000	