



PRODUCT DESCRIPTION

ProteinTech Group, Inc
2201 W. Campbell Park Dr.
Chicago, IL 60612
Toll Free: 1-888-4PTGLAB
In Illinois: 1312-455-8498
Fax: 1-312-455-8408
Web: www.ptglab.com

FOR IN VITRO RESEARCH USE ONLY
NOT FOR USE IN HUMANS OR ANIMALS

Antibody to MBP

Catalog Number: 66003-1-Ig

Isotype: IgG2a

Source: Mouse

Product Information

Description: Purified Mouse Anti-MBP monoclonal antibody

Immunogen: Recombinant Protein

Clonality: Monoclonal

Purify: Protein A Affinity Purified from mouse ascites fluid

Form: Liquid

Storage Condition: -20°C, Avoid freeze / thaw cycles.

Storage Buffer: PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3.

Concentration: 0.48mg/ml

Size: 100ul

Applications


Tested Applications: WB; IP; ELISA

Recommended Titers: ELISA 1:10,000~100,000

Western blotting 1:5,000~10,000

Specificity: Reacts with Escherichia coli.

Positive control: MBP-tagged protein



MBP-tagged fusion protein was subjected to SDS-PAGE followed by Western Blot with 66003-1-Ig at dilution of 1:8000.