

For Research Use Only

TUBGCP5 Polyclonal antibody

Catalog Number: 14620-1-AP **1 Publications**



Basic Information

Catalog Number:

14620-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6248

GenBank Accession Number:

BC071560

GeneID (NCBI):

114791

UNIPROT ID:

Q96RT8

Full Name:

tubulin, gamma complex associated protein 5

Calculated MW:

118 kDa

Observed MW:

118-120 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:400

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

mouse

Positive Controls:

WB : mouse brain tissue,

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Zhigang Li	31209365	Cell Death Differ	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse brain tissue were subjected to SDS PAGE followed by western blot with 14620-1-AP (TUBGCP5 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.