

For Research Use Only

MRPS10 Polyclonal antibody

Catalog Number: 16030-1-AP



Basic Information

Catalog Number:

16030-1-AP

Size:

200 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8897

GenBank Accession Number:

BC012560

GeneID (NCBI):

55173

UNIPROT ID:

P82664

Full Name:

mitochondrial ribosomal protein S10

Calculated MW:

201 aa, 23 kDa

Observed MW:

23 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:3000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HeLa cells, human liver tissue, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 16030-1-AP (MRPS10 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.