

For Research Use Only

SDR42E1 Polyclonal antibody

Catalog Number: 17347-1-AP



Basic Information

Catalog Number:

17347-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11391

GenBank Accession Number:

BC019670

GeneID (NCBI):

93517

UNIPROT ID:

Q8WUS8

Full Name:

short chain dehydrogenase/reductase
family 42E, member 1

Calculated MW:

383 aa, 43 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human bladder tissue slide using 17347-1-AP (SDR42E1 antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human bladder tissue slide using 17347-1-AP (SDR42E1 antibody) at dilution of 1:200 (under 40x lens).