

For Research Use Only

DDX49 Polyclonal antibody

Catalog Number: 20523-1-AP



Basic Information

Catalog Number:

20523-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14373

GenBank Accession Number:

BC002674

GeneID (NCBI):

54555

UNIPROT ID:

Q9Y6V7

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 49

Calculated MW:

483 aa, 54 kDa

Observed MW:

54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, Y79 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 20523-1-AP (DDX49 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.