

For Research Use Only

Advillin Polyclonal antibody

Catalog Number: 20956-1-AP **1 Publications**



Basic Information

Catalog Number: 20956-1-AP	GenBank Accession Number: BC111730	Purification Method: Antigen affinity purification
Size: 500 µg/ml	GeneID (NCBI): 10677	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: O75366	
Isotype: IgG	Full Name: advillin	
Immunogen Catalog Number: AG15114	Calculated MW: 819 aa, 92 kDa	
	Observed MW: 92 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse skin tissue, mouse small intestine tissue
Cited Applications: IF	
Species Specificity: human, mouse	
Cited Species: mouse	

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Charlotte Decourt	37068088	PLoS Biol	IF


Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 20956-1-AP (Advillin antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.