

For Research Use Only

Advillin Polyclonal antibody

Catalog Number: 20956-1-AP

1 Publications



Basic Information

Catalog Number:

20956-1-AP

Size:

500 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15114

GenBank Accession Number:

BC111730

GeneID (NCBI):

10677

UNIPROT ID:

O75366

Full Name:

advillin

Calculated MW:

819 aa, 92 kDa

Observed MW:

92 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

IF

Species Specificity:

human, mouse

Cited Species:

mouse

Positive Controls:

WB : mouse skin tissue, mouse small intestine tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Charlotte Decourt	37068088	PLoS Biol	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse skin tissue were subjected to SDS PAGE followed by western blot with 20956-1-AP (Advillin antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.