

For Research Use Only

CENPF Polyclonal antibody

Catalog Number: 20982-1-AP **3 Publications**



Basic Information

Catalog Number:

20982-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_016343

GeneID (NCBI):

1063

UNIPROT ID:

P49454

Full Name:

centromere protein F, 350/400ka
(mitosin)

Calculated MW:

368 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Cited Applications:

IF, IHC, WB

Species Specificity:

human

Cited Species:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Xuan Yang	31066149	Cancer Sci	WB
Min Wu	30655516	Nat Commun	IF
Dan Yang	36644760	Front Genet	IHC,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data