

For Research Use Only

C1orf92 Polyclonal antibody

Catalog Number: 21031-1-AP



Basic Information

Catalog Number: 21031-1-AP	GenBank Accession Number: BC033790	Purification Method: Antigen affinity purification
Size: 300 µg/ml	GeneID (NCBI): 149499	Recommended Dilutions: WB 1:200-1:1000 IF 1:20-1:200
Source: Rabbit	UNIPROT ID: Q8N4P6	
Isotype: IgG	Full Name: chromosome 1 open reading frame 92	
Immunogen Catalog Number: AG15344	Calculated MW: 559 aa, 62 kDa	
	Observed MW: 62 kDa	

Applications

Tested Applications: IF/ICC, WB, ELISA	Positive Controls:
Species Specificity: human	WB : HeLa cells, IF : HeLa cells,

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 21031-1-AP (C1orf92 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of HeLa cells using 21031-1-AP (C1orf92 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).