

For Research Use Only

# TRA2B Polyclonal antibody

Catalog Number: 23832-1-AP **1 Publications**



## Basic Information

<b>Catalog Number:</b> 23832-1-AP	<b>GenBank Accession Number:</b> BC005898	<b>Purification Method:</b> Antigen affinity purification
<b>Size:</b> 230 µg/ml	<b>GeneID (NCBI):</b> 6434	
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P62995	
<b>Isotype:</b> IgG	<b>Full Name:</b> transformer 2 beta homolog (Drosophila)	
<b>Immunogen Catalog Number:</b> AG20683	<b>Calculated MW:</b> 34 kDa	
	<b>Observed MW:</b> 34 kDa	

## Applications

**Tested Applications:**  
ELISA

**Cited Applications:**  
WB

**Species Specificity:**  
human, mouse

**Cited Species:**  
human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Shu-Qi Li	32405436	J Adv Res	WB

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data