

For Research Use Only

MAGEF1 Polyclonal antibody

Catalog Number: 24105-1-AP



Basic Information

Catalog Number:	24105-1-AP	GenBank Accession Number:	BC010056	Purification Method:	Antigen affinity purification
Size:	700 µg/ml	GeneID (NCBI):	64110	Recommended Dilutions:	IHC 1:20-1:200
Source:	Rabbit	UNIPROT ID:	Q9HAY2		
Isotype:	IgG	Full Name:	melanoma antigen family F, 1		
Immunogen Catalog Number:	AG21147	Calculated MW:	35 kDa		

Applications

Tested Applications:	IHC, ELISA	Positive Controls:	IHC : human testis tissue,
Species Specificity:	human		

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 24105-1-AP (MAGEF1 Antibody) at dilution of 1:162 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 24105-1-AP (MAGEF1 Antibody) at dilution of 1:162 (under 40x lens).