

For Research Use Only

# MAGEF1 Polyclonal antibody

Catalog Number: 24105-1-AP



## Basic Information

Catalog Number: 24105-1-AP	GenBank Accession Number: BC010056	Purification Method: Antigen affinity purification
Size: 700 µg/ml	GeneID (NCBI): 64110	Recommended Dilutions: IHC 1:20-1:200
Source: Rabbit	UNIPROT ID: Q9HAY2	
Isotype: IgG	Full Name: melanoma antigen family F, 1	
Immunogen Catalog Number: AG21147	Calculated MW: 35 kDa	

## Applications

Tested Applications:  
IHC, ELISA  
Species Specificity:  
human

Positive Controls:  
IHC : human testis tissue,

**Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (\*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0**

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.  
Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 24105-1-AP (MAGEF1 Antibody) at dilution of 1:162 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human testis tissue slide using 24105-1-AP (MAGEF1 Antibody) at dilution of 1:162 (under 40x lens).