

For Research Use Only

ABTB2 Polyclonal antibody

Catalog Number: 24323-1-AP



Basic Information

Catalog Number:

24323-1-AP

Size:

750 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19284

GenBank Accession Number:

BC093968

GeneID (NCBI):

25841

UNIPROT ID:

Q8N961

Full Name:

ankyrin repeat and BTB (POZ) domain containing 2

Calculated MW:

839 aa, 93 kDa

Observed MW:

114 kDa, 93 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HepG2 cells, K-562 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



HepG2 cells were subjected to SDS PAGE followed by western blot with 24323-1-AP (ABTB2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.