

For Research Use Only

MRPL47 Polyclonal antibody

Catalog Number: 24728-1-AP



Basic Information

Catalog Number: 24728-1-AP	GenBank Accession Number: BC032522	Purification Method: Antigen affinity purification
Size: 1000 µg/ml	GeneID (NCBI): 57129	Recommended Dilutions: WB 1:200-1:1000
Source: Rabbit	UNIPROT ID: Q9HD33	
Isotype: IgG	Full Name: mitochondrial ribosomal protein L47	
Immunogen Catalog Number: AG11744	Calculated MW: 230 aa, 27 kDa	
	Observed MW: 29 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse testis tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 24728-1-AP (MRPL47 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.