

For Research Use Only

ABHD1 Polyclonal antibody

Catalog Number: 26091-1-AP



Basic Information

Catalog Number: 26091-1-AP	GenBank Accession Number: BC039576	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 84696	Recommended Dilutions: WB 1:200-1:500
Source: Rabbit	UNIPROT ID: Q965E0	
Isotype: IgG	Full Name: abhydrolase domain containing 1	
Immunogen Catalog Number: AG23138	Calculated MW: 405 aa, 45 kDa	
	Observed MW: 48 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse heart tissue,
Species Specificity: human, mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 26091-1-AP (ABHD1 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.