

For Research Use Only

ABHD1 Polyclonal antibody

Catalog Number: 26091-1-AP



Basic Information

Catalog Number:

26091-1-AP

Size:

450 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23138

GenBank Accession Number:

BC039576

GeneID (NCBI):

84696

UNIPROT ID:

Q96SE0

Full Name:

abhydrolase domain containing 1

Calculated MW:

405 aa, 45 kDa

Observed MW:

48 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:200-1:500

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse heart tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 26091-1-AP (ABHD1 Antibody) at dilution of 1:300 incubated at 4 degree celsius over night.