

For Research Use Only

CCNYL1 Polyclonal antibody

Catalog Number: 26970-1-AP



Basic Information

Catalog Number:

26970-1-AP

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25583

GenBank Accession Number:

BC067253

GeneID (NCBI):

151195

UNIPROT ID:

Q8N7R7

Full Name:

cyclin Y-like 1

Calculated MW:

238 aa, 28 kDa

Observed MW:

45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : fetal human brain tissue, THP-1 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 26970-1-AP (CCNYL1 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.