

For Research Use Only

# ABI3BP Polyclonal antibody

Catalog Number: 27727-1-AP



## Basic Information

Catalog Number:

27727-1-AP

Size:

800 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26841

GenBank Accession Number:

BC030221

GeneID (NCBI):

25890

UNIPROT ID:

Q7Z7G0

Full Name:

ABI family, member 3 (NESH) binding protein

Calculated MW:

119 kDa

Observed MW:

54 kDa, 72 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, Mouse

Positive Controls:

WB: C2C12 cell, mouse thymus tissue, C2C12 cells, mouse testis tissue

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27727-1-AP (ABI3BP antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.