

For Research Use Only

FAM55B Polyclonal antibody

Catalog Number: 27877-1-AP



Basic Information

Catalog Number:

27877-1-AP

Size:

900 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27351

GenBank Accession Number:

BC156724

GeneID (NCBI):

120406

UNIPROT ID:

Q96DL1

Full Name:

family with sequence similarity 55, member B

Calculated MW:

559 aa, 65 kDa

Observed MW:

60-65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:1000-1:8000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HepG2 cells, MCF-7 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 27877-1-AP (FAM55B antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.