

For Research Use Only

ZNHIT6 Polyclonal antibody

Catalog Number: 28019-1-AP



Basic Information

Catalog Number: 28019-1-AP	GenBank Accession Number: BC026236	Purification Method: Antigen affinity purification
Size: 900 µg/ml	GeneID (NCBI): 54680	Recommended Dilutions: WB 1:1000-1:4000
Source: Rabbit	UNIPROT ID: Q9NWK9	
Isotype: IgG	Full Name: zinc finger, HIT type 6	
Immunogen Catalog Number: AG27792	Calculated MW: 470 aa, 54 kDa	
	Observed MW: 50-55 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : mouse lung tissue,
Species Specificity: Human, Mouse	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 28019-1-AP (ZNHIT6 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.