

For Research Use Only

ST6GAL2 Polyclonal antibody

Catalog Number: 28367-1-AP **1 Publications**



Basic Information

Catalog Number:

28367-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26357

GenBank Accession Number:

BC008680

GeneID (NCBI):

84620

UNIPROT ID:

Q96JFO

Full Name:

ST6 beta-galactosamide alpha-2,6-sialyltransferase 2

Calculated MW:

529 aa, 60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:800

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse

Cited Species:

human

Positive Controls:

WB : HL-60 cells, RAW 264.7 cells, K-562 cells

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Ruijia Liu	38318244	J Inflamm Res	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28367-1-AP (ST6GAL2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.