

For Research Use Only

# ST6GAL2 Polyclonal antibody

Catalog Number: 28367-1-AP

1 Publications



## Basic Information

Catalog Number:

28367-1-AP

Size:

850 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26357

GenBank Accession Number:

BC008680

GeneID (NCBI):

84620

UNIPROT ID:

Q96JFO

Full Name:

ST6 beta-galactosamide alpha-2,6-sialyltransferase 2

Calculated MW:

529 aa, 60 kDa

Observed MW:

60 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:300-1:800

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse

Cited Species:

human

Positive Controls:

WB : HL-60 cells, RAW 264.7 cells, K-562 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Ruijia Liu	38318244	J Inflamm Res	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28367-1-AP (ST6GAL2 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.