

For Research Use Only

GABRB1 Polyclonal antibody

Catalog Number: 28425-1-AP **1 Publications**



Basic Information

Catalog Number:

28425-1-AP

Size:

700 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29243

GenBank Accession Number:

BC022449

GeneID (NCBI):

2560

UNIPROT ID:

P18505

Full Name:

gamma-aminobutyric acid (GABA) A receptor, beta 1

Calculated MW:

474 aa, 54 kDa

Observed MW:

50-54 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

Human, Mouse, Rat

Cited Species:

human

Positive Controls:

WB : mouse brain tissue, rat brain tissue

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Min Wang	38292191	Genes Dis	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28425-1-AP (GABRB1 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.