

For Research Use Only

TCTN1 Polyclonal antibody

Catalog Number: 28630-1-AP



Basic Information

Catalog Number: 28630-1-AP	GenBank Accession Number: BC062611	Purification Method: Antigen affinity purification
Size: 450 µg/ml	GeneID (NCBI): 79600	Recommended Dilutions: WB 1:500-1:1000
Source: Rabbit	UNIPROT ID: Q2MV58	
Isotype: IgG	Full Name: tectonic family member 1	
Immunogen Catalog Number: AG30064	Calculated MW: 64 kDa	
	Observed MW: 60-75 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : DU 145 cells, mouse brain tissue, HCT 116 cells
Species Specificity: Human, mouse, rat	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 28630-1-AP (TCTN1 antibody) at dilution of 1:800 incubated at room temperature for 1.5 hours.