

For Research Use Only

KIAA0284 Monoclonal antibody

Catalog Number: 66749-1-Ig



Basic Information

Catalog Number:

66749-1-Ig

Size:

1800 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26852

GenBank Accession Number:

BC112928

GeneID (NCBI):

283638

UNIPROT ID:

Q9Y4F5

Full Name:

KIAA0284

Calculated MW:

1589 aa, 172 kDa

Observed MW:

172 kDa

Purification Method:

Protein A purification

CloneNo.:

1D4B8

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB : HeLa cells, HEK-293 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66749-1-Ig (KIAA0284 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.