

ARF5 Monoclonal antibody

Catalog Number: 67279-1-Ig

Basic Information

Catalog Number:

67279-1-Ig

Size:

2000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7849

GenBank Accession Number:

BC003043

GeneID (NCBI):

381

UNIPROT ID:

P84085

Full Name:

ADP-ribosylation factor 5

Calculated MW:

21 kDa

Observed MW:

18-20 kDa

Purification Method:

Protein G purification

CloneNo.:

2B3B1

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat, pig

Positive Controls:

WB : A549 cells, pig brain tissue, rat brain tissue, mouse brain tissue, HeLa cells, HepG2 cells, Jurkat cells

Background Information

Storage

Storage:


Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67279-1-Ig (ARF5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.