

For Research Use Only

BRD4 Monoclonal antibody

Catalog Number: 67374-1-Ig **1 Publications**



Basic Information

Catalog Number: 67374-1-Ig	GenBank Accession Number: BC035266	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 23476	CloneNo.: 1F11A6
Source: Mouse	UNIPROT ID: O60885	
Isotype: IgG1	Full Name: bromodomain containing 4	
Immunogen Catalog Number: AG27545	Calculated MW: 1362 aa, 152 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications:
ELISA

Cited Applications:
IF, WB

Species Specificity:
Human

Cited Species:
mouse

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Yanwen Lan	32596881	J Cell Mol Med	WB,IF

Storage

Storage:
Store at -20°C.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data