

For Research Use Only

# SALL4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83039-1-PBS

## Basic Information

Catalog Number:

83039-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16076

GenBank Accession Number:

BC111714

GeneID (NCBI):

57167

UNIPROT ID:

Q9UJQ4

Full Name:

sal-like 4 (Drosophila)

Calculated MW:

1053 aa, 112 kDa

Purification Method:

Protein A purification

CloneNo.:

23028487

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00137-1 measured by sandwich ELISA. Capture antibody: 83039-2-PBS. Detection antibody: 83039-1-PBS. Range: 0.156-10 ng/mL