

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Cd2ap Fusion Protein



## Basic Information

**Catalog Number:**  
Ag0622

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
mouse

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Validated Application:**  
Blocking peptide

**Peptide Sequence:**  
MVDYIVEYDYDAVHDELTIKRVGEIIRNVKKLQEEGW  
LEGELNGRRGMFPDNFVKEIKRETEPKDDNLPKREKQ  
GNVASLVQRISTYGLPAGGIQPHPQTKAIKKTKKRQ  
CKVLFDYSPQNEDELELIVGDVIDVIEVEEGWWSGT  
LNNKLGFLPSNFVKELESTEDGETHNAQESEVPLTGP  
TSPSPGNGSEAPGSVAQPKKIRGIGFDIFKEGSV  
KLRTRTSSSETEEKTEKPLILQPLGSRQNVVETKPDV  
DGKIKAKEYCRTLFPYGTNEDELTFREGEIHLISKETG  
EAGWWWKELNGKEGVFPDNFAVQISELDKDFPKPKK  
PPPPAKGPAPKPDLSAAEKKAFPLKAEKDEKSLLEQK  
PSKPAAPQVPPKPTAPTASNLLRSPGAVYVYKRPKPK  
VPPPPAAKINGEVSISSKIDTEPVSKPKLDPQLPVR  
PKSVDLDAFVARNKSKETDDVNFDDIASSENLLHTAN  
RPKMPGRRRLPGRFNGGHSPTQSPEKTLKLPKEDDSGNL  
KPLEFKKDA SYSSKPSLSTPSSASKVNTAAFLTPELKA  
KA EADDGKKNSVDELRAQIIELLCIVDALKKDHGKELE  
KLRKELEEEKAMRSNLEVEIAKLLKAVLLS  
(1-637 aa encoded by BC157901)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.