

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TORC1 Fusion Protein



## Basic Information

**Catalog Number:**

Ag0693

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MATSNNPRKFSEKIALHNQKQAEETA AFEEVMKDLSL  
TRAA RLQ LQK SQYLQ L G P S R G Q Y Y G G S L P N V N Q I G S  
G T M D L P F Q T P F Q S S G L D T S R T T R H H G L V D R V Y R E R G R  
L G S P H R R P L S V D K H G R Q A D S C P Y G T M Y L S P P A D T S W  
R R T N S D S A L H Q S T M T P T Q P E S F S S G S Q D V H Q K R V L L L  
T V P G M E E T T S E A D K N L S K Q A W D T K K T G S R P K S C E V P  
G I N I F P S A D Q E N T T A L I P A T H N T G G S L P D L T N I H F P S P L  
P T P L D P E E P T F P A L S S S S T G N L A A N L T H L G I G G A G Q G  
M S T P

(1-300 aa encoded by BC028050)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.