

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Cd2ap Fusion Protein



Basic Information

Catalog Number:
Ag21170

Size:
50 µg

Form:
Available lyophilized

Species:
mouse

Expression Source:
e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:
MVDYIVEYDYDAVHDELTI RVGEIIRNVKKLQEEGW
LEGELNGRRGMFPDNFVKEIKRETEPKDDNLPKREERQ
GNVASLVQRISTYGLPAGGIQPHPQTKAIKKTKKRQ
CKVLFDYSPQNEDELELIVGDVIDVIEVEEGWWSGT
LNNKLG LFP SNFVKELESTEDGETHNAQESEVPLTGP
TSP LPSPGNGSEPA PGSVAQPKKIRGIGFDIFKEGSV
KLRTRTSSSETEEEKTEKPLILQPLGSR TQNVEVTKPDV
DGKIKAKEYCRTLFPYTG TNEDELTFREGEIHLISKETG
EAGWWWGELNGKEGVFPDNFAVQISELDKDFPKPKK
PPPPAKGPAPKPDLSAAEKKAFPLKAEKDEKSLLEQK
PSKPAAPQVPPKPTAPT KASNLLRSPGAVYPKRPEKP
VPPPPAAKINGEVSISSKIDTEPVSKPKLDPEQLPVR
PKSVDLDAFVARN SKETDDVNFDDIASSENLLHTAN
RPKMPGRRLPGRFNGGHSPTQSPEKTLKLPKEDDSGNL
KPLEFKKDA SYSSKPSLSTPSSASKVNTAAFLTPELKA
KA EADDGKKNSVDELRAQI IELLCIVDALKKDHGKELE
KLRKELEEEKAMRSNLEVEIAKLKAVLLS
(1-637 aa encoded by BC157901)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.