

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

DMAP1 Fusion Protein



Basic Information

Catalog Number:

Ag28906

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MATGADVVDILELGGPEGDAASGTISKKDIINPDKKKS
KKSETLTFKRPEGMHREYVYALLYSDKIDAPLLPSDT
GQGYRTVKAKLGSKKVRPWKWPFTNPARKDGAMF
FHWRRAAEEGKDYPFARFNKTVQVPVYSEQEYQLYL
HDDAWTKAETDHLFDLSRRFDLRFVVIHDRYDHQQF
KKRSVDELKERYHYHICAKLANVRAVPGTDLKIPVFDAG
HERRRKEQLERLYNRTPEQVAEEYLLQELRKIEARKKE
REKRSQDLQKLITAAADTTAEQRRTERKAPKKLPQKKE
AEKPAVPETAGIKFPDFKSAAGVTLRSQRMKLPSSVGQ
KKIKALEQMELLELGVLSPTTEELVHMFNELRSDLVLL
YELKQACANCEYELQMLRHRHEALARAGVLGGPATP
ASGPGPASAEPAVTEPGLGDPKDTIIDVVGAPLTPN
SRKRRESASSSSVKKAKKP
(1-467 aa encoded by BC002855)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.