

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# 143cn Fusion Protein



## Basic Information

**Catalog Number:**

Ag3909

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

**Expression Source:**

*e coli*-derived, PET28a, with N-terminal 6<sup>°</sup>His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

MNRVKVAQDLVNERLYRTRPEALSFTSIVDTPEVVLA  
KANSLQISEKLYQEAWNKDKSNITIPSDTPEMLQAHIN  
ALQISNKLYQKDWNDTKQKGYDIRADAIEIKHAKASR  
EIASEYKYKEGYRKQLGHHMGFRTLQDDPKSVWAIHA  
AKIQSDREYKKAYEKSKGIHNTPLDMMMSIVQAKKCQV  
LVSDIDYRNYLHQWTCCLPDQNDVIQAKKAYDLQSDN  
LYKSDLEWLKIGWLPEGSVEVMRVKNAQNLNERLY  
RIKPEALKFTSIVDTPEVIQAKINAVQISEPLY  
(1-235 aa encoded by )

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.