

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# ATP2A3 Fusion Protein



## Basic Information

**Catalog Number:**

Ag4527

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

TRRMARKNAIVRSLPSVETLGCTSVICSDKTGTLTTNQ  
MSVCRMVVAEADAGSCLLHEFTISGTTYTPEGEVR  
QGDDQPVRCGQFDGLVELATICALCNDSDYNEAKG  
VYEKVGAEATETALTCLVEKMNVDLQALSRVERAG  
ACNTVIKQLMRKEFTLEFSRDRKSMVYCTPTRPHPT  
GQGSKMVFVKAPESVIERCSSVRVGSRTAPLTPSRE  
QILAKIRDWGSGLRCLALATRDAPPRKEDMELDD  
CSKFVQYETDLTFVGCVGMLDPPRPEVAACITRCYQ  
AGIRVVMITGDNKGTAVAICRRLGIFGDTEVAGKAY  
TGREFDLSPQQACRTARCFARVEPAHKSRIVE  
NLQSFNEITAMTGDGVNDAPALKKAEIGIAMGSGTA  
VAKSAAEMVLSDDNFASIVAAVEEGRAIYSNMKQFIR  
YLISSNVGEVVCIFLTAIILGLPEALIPVQLLWNLVTDG  
LPATALGFNPPDLIMEKLPRSPREALISGWLFRRYLAI  
GVVYGLATVAAAATWVYDAEGPHINFYQLRNFLKC  
SEDNPLFAGIDCEVFESRFPTTMALSVLVTIEMCNALN  
SVSENQSLLRMPWMPWLLVAVAMSMALHFLILLV  
PPLPLIFQVTPLSGRQWVVVLIQISLPVILLDEALKYLSR  
NHMHEEMSQK

(323-999 aa encoded by BC035729)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.