

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

FUCA1 Fusion Protein



Basic Information

Catalog Number:

Ag9612

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

EEWADLFQAAGAKYVVLTTKHHEGFTNWPSVSWN
WNSKDVGPHRDLV GELGTALRKRNI RYGLYHSLEWF
HPLYLLDKKNGFKTQHFVSAKTMPELYDLVNSYKPLI
WSDGEWECPDYWNSTNFLSWLYNDSPVKDEVVVN
DRWGQNCSCHHGGYNCEDKFKPQSLPDHKWEMC
TSIDKFSWGYRRDMALSDVTEESEIISELVQTVSLGGN
YLLNIGPTKDGLIVPIFQERLLAVGKWL SINGEAIYASK
PWRVQWEKNTTSVWYTSKGS AVYAI FLHWPENGL
NLESPITTTTKITMLGIQ GDLKWSTDPDKGLFISLPQL
PPSAVPAEFAWTIKLTGVK

(113-461 aa encoded by BC017338)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.