

For Research Use Only

# CoraLite® Plus 488-conjugated CNBP Monoclonal antibody



Catalog Number: CL488-67109

## Basic Information

Catalog Number:

CL488-67109

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG7114

GenBank Accession Number:

BC000288

GeneID (NCBI):

7555

UNIPROT ID:

P62633

Full Name:

CCHC-type zinc finger, nucleic acid binding protein

Calculated MW:

19 kDa

Observed MW:

19-25 kDa

Purification Method:

Protein A purification

CloneNo.:

1H1D12

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC (Intra)

Species Specificity:

Human, mouse, rat

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> HeLa cells were intracellularly stained with 0.4 ug Coralite® Plus 488 Anti-Human CNBP (CL488-67109, Clone:1H1D12) (red), or 0.4 ug Mouse IgG2a Isotype Control (CL488-66360-2, Clone: K11A1B2A2) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).