

For Research Use Only.
Not For Use In Diagnostics.

CoraLite647-conjugated AffiniPure F(ab')₂ Fragment Goat Anti-Mouse IgG (H+L)



Catalog Number:SA00014-10

Information

Catalog Number: SA00014-10	Reactivity: Mouse
Host: Goat	Physical State: Liquid
Applications: IF, FC	Conjugation: CoraLite®647

Recommended Dilutions

1:100 – 1:500

Preservative

0.01% Sodium Azide

Fluorophore

CL647 Fluorescent Dye. A_{max}=642nm, E_{max}=662nm

Stabilizer

5 mg/ml BSA

Absorbance Ratio

A₆₅₁/A₂₈₀ = 4.6

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 7.6.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data